

NIBAPLEX®



Off-the-shelf solutions
for viral vector
characterisation

MULTIPLEX AND SIMPLEX ASSAYS FOR VECTOR GENOME INTEGRITY & TITRE

NibaPlex® digital PCR (dPCR) assays have been developed for specific detection and quantification of a target panel of viral vector genome elements, empowering scientists in biopharma industry to address genome integrity and vector genome titre.



1 VECTOR GENOME INTEGRITY

Cutting-edge multiplex dPCR
assays for genome integrity
assessment

The key role of multiplex dPCR for the quantification of complete genomes in your viral vector product

Genome integrity results by multiplex dPCR show an excellent correlation of complete vector genomes to the results of biopotency, as presented in recent studies.

Our dedicated NibaPlex® multiplex dPCR assays have been developed for specific detection and quantification of complete and incomplete genomes, providing in-depth information on the genome integrity of viral vector products for gene therapy.

Our portfolio of NibaPlex® assays can provide detailed information on the quantity of complete genomes and the distribution of incomplete populations in AAV samples.



BENEFITS


High throughput method to help you improve biopotency of viral vector products.

NIBAPLEX® 4-PLEX ASSAYS

First-of-its-kind assay with
4 target regions
simultaneously targeted
in one reaction mix

The NibaPlex® 4-plex assays are designed to simultaneously target 4 regions on the viral vector genome to obtain in-depth information on genome integrity. This enables the user to collect quantitative data of complete and incomplete vector genomes present in viral vector capsids and to determine the expected biopotency of viral vector products.

NibaPlex® 4-plex off-the-shelf offerings comprise various combinations of enhancer, promoter, and poly-A tail options with the model GFP gene.



| NibaPlex® 4-plex assay | 1 Enhancer | 2 Promoter | 3 GOI | 4 poly-A tail |
|------------------------|------------|------------|-------|---------------|
| 4.01 | CMV | CMV | GFP | SV40 |
| 4.02 | CMV | CMV | GFP | hGH |
| 4.03 | CMV | CMV | GFP | bGH |



CUSTOM DESIGN & DEVELOPMENT OF 4-PLEX ASSAYS


Niba Labs offers custom design of 4-plex assays for your specific regions & GOI upon request.

NIBAPLEX® DUPLEX ASSAYS

First-stage assessment of viral vector genome integrity

The NibaPlex® duplex assays are designed to simultaneously target the beginning and end of the viral vector genome to obtain high-level information on genome integrity based on quantitative data of complete and incomplete vector genomes present in viral vector capsids. This enables the user to quickly assess the general efficacy of the viral vector production process.

NibaPlex® duplex off-the-shelf offerings comprise various combinations of enhancer and poly-A tail options.



| NibaPlex® duplex assay | 1 Enhancer | 2 Promoter | 3 GOI | 4 poly-A tail |
|------------------------|------------|------------|-------|---------------|
| 2.01 | CMV | - | - | SV40 |
| 2.02 | CMV | - | - | hGH |
| 2.03 | CMV | - | - | bGH |



CUSTOM DESIGN & DEVELOPMENT OF DUPLEX ASSAYS


Niba Labs offers custom design of duplex assays for your specific vector genome sequence.

2 VECTOR GENOME TITRE

Simplex dPCR assays for quantification of viral vector genomes




The NibaPlex® simplex assays are highly-optimized to target a single element of the viral vector genome, enabling the user to determine the titer of the viral genomes with great accuracy and precision. Our broad portfolio includes assays targeting the most common viral vector genome elements, rendering them appropriate for the quantification of commercial viral vector products.



| NibaPlex® simplex assay | 1 Enhancer | 2 Promoter | 3 GOI | 4 poly-A tail |
|-------------------------|------------|------------|-------|---------------|
| 1.01 | CMV | - | - | - |
| 1.02 | - | CMV | - | - |
| 1.03 | - | - | GFP | - |
| 1.04 | - | - | - | SV40 |
| 1.05 | - | - | - | hGH |
| 1.06 | - | - | - | bGH |

3 PROCESS-RELATED IMPURITIES

Simplex and multiplex dPCR assays for analysis of impurities in viral vector products



The NibaPlex® Kan^R assays are specifically designed to enable the accurate determination of the presence, quantity, and fragmentation of residual Kan^R gene in all stages of the viral vector production process.

The NibaPlex® Kan^R multiplex dPCR assays are based on an innovative approach for assessment of complete and incomplete populations of the plasmid-derived kanamycin resistance gene (nptII). The duplex assay targeting the start and end of the Kan^R gene provides high-level information on the complete and incomplete gene sequences, with the triplex assay enabling further assessment of their fragmentation. Our NibaPlex® Kan^R simplex assay targeting the middle of the Kan^R gene enables its general titration.

| NibaPlex® Kan ^R assay | Assay | Use | Kan ^R Target 1 | Kan ^R Target 2 | Kan ^R Target 3 |
|----------------------------------|---------|---|---------------------------|---------------------------|---------------------------|
| 1.07 | simplex | Quantification of Kan ^R gene | - | Middle of gene | - |
| 2.04 | duplex | High-level determination of Kan ^R gene integrity | Start of gene | - | End of gene |
| 3.01 | triplex | In-depth determination of Kan ^R gene integrity | Start of gene | Middle of gene | End of gene |

OFF-THE-SHELF NIBAPLEX[®] ASSAYS

Exceptionally user friendly:
optimised assays, ready-to-use &
all reagents in one tube

Assays are provided as single tube containing ready-to-use, 20x-concentrated reference assay

| Cat. No. / ID | Item | Type of assay | Target region | Dye | Number of reactions |
|---------------|---------------|----------------------------|--|--------------------------------|--|
| 4-PLEX | | | | | |
| 400011 | NibaPlex 4.01 | Vector Genome Integrity | CMV enhancer CMV promoter GFP SV40 poly-A | FAM Cy5 Texas Red HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 400015 | NibaPlex 4.01 | Vector Genome Integrity | CMV enhancer CMV promoter GFP SV40 poly-A | FAM Cy5 Texas Red HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 400021 | NibaPlex 4.02 | Vector Genome Integrity | CMV enhancer CMV promoter GFP hGH poly-A | FAM Cy5 Texas Red HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 400025 | NibaPlex 4.02 | Vector Genome Integrity | CMV enhancer CMV promoter GFP hGH poly-A | FAM Cy5 Texas Red HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 400031 | NibaPlex 4.03 | Vector Genome Integrity | CMV enhancer CMV promoter GFP bGH poly-A | FAM Cy5 Texas Red HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 400035 | NibaPlex 4.03 | Vector Genome Integrity | CMV enhancer CMV promoter GFP bGH poly-A | FAM Cy5 Texas Red HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| TRIPLEX | | | | | |
| 300011 | NibaPlex 3.01 | Process-Related Impurities | Kanamycin (three regions) | FAM Cy5 HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 300015 | NibaPlex 3.01 | Process-Related Impurities | Kanamycin (three regions) | FAM Cy5 HEX | sufficient for 5x96 dPCR reactions of 12 µl each |

OFF-THE-SHELF NIBAPLEX[®] ASSAYS

Exceptionally user friendly:
optimised assays, ready-to-use &
all reagents in one tube

Assays are provided as single tube containing ready-to-use, 20x-concentrated reference assay

| Cat. No. / ID | Item | Type of assay | Target region | Dye | Number of reactions |
|---------------|---------------|----------------------------|---------------------------------------|------------|--|
| DUPLEX | | | | | |
| 200011 | NibaPlex 2.01 | Vector Genome Integrity | CMV enhancer SV40 poly-A | FAM HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 200015 | NibaPlex 2.01 | Vector Genome Integrity | CMV enhancer SV40 poly-A | FAM HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 200021 | NibaPlex 2.02 | Vector Genome Integrity | CMV enhancer hGH poly-A | FAM HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 200025 | NibaPlex 2.02 | Vector Genome Integrity | CMV enhancer hGH poly-A | FAM HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 200031 | NibaPlex 2.03 | Vector Genome Integrity | CMV enhancer bGH poly-A | FAM HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 200035 | NibaPlex 2.03 | Vector Genome Integrity | CMV enhancer bGH poly-A | FAM HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 200041 | NibaPlex 2.04 | Process-Related Impurities | Kanamycin (beginning and end of gene) | FAM HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 200045 | NibaPlex 2.04 | Process-Related Impurities | Kanamycin (beginning and end of gene) | FAM HEX | sufficient for 5x96 dPCR reactions of 12 µl each |

OFF-THE-SHELF NIBAPLEX[®] ASSAYS

Exceptionally user friendly:
optimised assays, ready-to-use &
all reagents in one tube

Assays are provided as single tube containing ready-to-use, 20x-concentrated reference assay

| Cat. No. / ID | Item | Type of assay | Target region | Dye | Number of reactions |
|---------------|---------------|----------------------------|---------------|-----------|--|
| SIMPLEX | | | | | |
| 100011 | NibaPlex 1.01 | Vector Genome Titer | CMV enhancer | FAM | sufficient for 1x96 dPCR reactions of 12 µl each |
| 100015 | NibaPlex 1.01 | Vector Genome Titer | CMV enhancer | FAM | sufficient for 5x96 dPCR reactions of 12 µl each |
| 100021 | NibaPlex 1.02 | Vector Genome Titer | CMV promoter | Cy5 | sufficient for 1x96 dPCR reactions of 12 µl each |
| 100025 | NibaPlex 1.02 | Vector Genome Titer | CMV promoter | Cy5 | sufficient for 5x96 dPCR reactions of 12 µl each |
| 100031 | NibaPlex 1.03 | Vector Genome Titer | GFP | Texas Red | sufficient for 1x96 dPCR reactions of 12 µl each |
| 100035 | NibaPlex 1.03 | Vector Genome Titer | GFP | Texas Red | sufficient for 5x96 dPCR reactions of 12 µl each |
| 100041 | NibaPlex 1.04 | Vector Genome Titer | SV40 poly-A | HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 100045 | NibaPlex 1.04 | Vector Genome Titer | SV40 poly-A | HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 100051 | NibaPlex 1.05 | Vector Genome Titer | hGH poly-A | HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 100055 | NibaPlex 1.05 | Vector Genome Titer | hGH poly-A | HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 100061 | NibaPlex 1.06 | Vector Genome Titer | bGH poly-A | HEX | sufficient for 1x96 dPCR reactions of 12 µl each |
| 100065 | NibaPlex 1.06 | Vector Genome Titer | bGH poly-A | HEX | sufficient for 5x96 dPCR reactions of 12 µl each |
| 100071 | NibaPlex 1.07 | Process-Related Impurities | Kanamycin | Cy5 | sufficient for 1x96 dPCR reactions of 12 µl each |
| 100075 | NibaPlex 1.07 | Process-Related Impurities | Kanamycin | Cy5 | sufficient for 5x96 dPCR reactions of 12 µl each |